



Product Manual

TaqGen® Blood Genomic DNA Kit 20 preps

Cat no: MSB120

The TaqGen® Blood Genomic DNA Kit is suitable for Genomic DNA Extraction from whole Blood, body fluids, buccal swabs, buffy coat and avian blood.

For research use only.

Not for use in diagnostic procedures.

Description

The TaqGen® Blood Genomic DNA kit, a spin column based kit, provide an easy, rapid and efficient purification of high quality genomic DNA from whole blood and body fluid. The purified DNA can be used directly in a variety of downstream applications, including PCR, RealTime PCR, southern blotting and restriction enzyme digestion. The kit eliminates the need for expensive resin, toxic phenol-chloroform extractions, or time-consuming alcohol precipitation. The standard procedure takes less than 20 minutes following cell lysis and yield purified DNA greater than 30 Kb in size. This kit allows for the single or multiple simultaneous processing of samples.

Product Components and Storage Conditions

| Product | Storage condition | Total Vol | Vol/Reaction For 200 µL sample |
|-----------------|-------------------|-----------|--------------------------------|
| Proteinase K | 2°-8°C | 400µL | 20 µL |
| Lysis Buffer | RT | 20 ml | 400 µL |
| Wash Buffer WB1 | RT | 8 ml conc | 500 µL |
| Wash Buffer WB2 | RT | 4 ml conc | 500 µL |
| Elution Buffer | RT | 20 ml | 100 µL |
| Spin Columns | RT | 20 pc | 1 |
| Collection Tube | RT | 20 pc | 1 |

Reconstitution of buffers/reagents:

| Buffer/Enzyme name | Procedure |
|--------------------------------|------------------------------|
| Proteinase K Lyophilized (5mg) | Add 200 µL Proteinase buffer |
| Wash Buffer WB1 (8 ml conc) | Add 6 ml 100% Ethanol |
| Wash Buffer WB2 (4 ml conc) | Add 10 ml 100% Ethanol |

Storage Conditions: Store the TaqGen® Blood Genomic DNA Kit at 15°-30°C except the proteinase K after reconstitute.

Safety Information: The TaqGen® Blood Genomic DNA kit contains guanidine hydrochloride and proteinase K, should be considered harmful and irritants.

Caution: Handle buffer with care. Bleach reacts with guanidine hydrochloride and should not be added to any sample waste.

Before You Begin

To minimize DNA degradation, avoid repeated freeze/thaw cycles of the samples and perform extractions from fresh material or material that has been immediately frozen and stored at -20 °C or -70 °C. Buffers For sample volume adjustment: PBS:137 mM NaCl, 2.7 mM KCl, 10 mM Na₂HPO₄, 2 mM KH₂PO₄, pH 7.4.

TE buffer: 10 mM Tris-HCl, pH 8.0, 1 mM EDTA.

Protocol for Whole Blood Genomic DNA Purification

1. Lyse Cells

- a. Add **20 μ L of Proteinase K Solution** to **200 μ L of whole blood**, mix by vortexing. Add **400 μ L of Lysis Buffer**, mix thoroughly by vortexing or pipetting to obtain a uniform suspension

Note. If using less than 200 μ L of blood, adjust sample volume to 200 μ L with 1X PBS or TE buffer (not provided).

- b. Incubate at 56°C for 10 min in a water bath. While vortexing occasionally or use a shaking water bath, rocking platform or thermomixer until the cells are completely lysed.

2. Adjust DNA binding conditions

- a. Add **200 μ L of ethanol** (96-100%).
Note: Adjust the volume of ethanol required based on the amount of starting material. Vortex to mix thoroughly

3. Bind DNA

- a. For each sample, place one **Spin Column** into a Collection Tube. Apply the sample to the column. Centrifuge for **1 min** at **11,000 x g**. Discard the filtrate and reuse the collection tube.

4. Wash and dry silica membrane

- a. Add **500 μ L Wash Buffer WB1** (with Ethanol Added). Centrifuge for **1 min** at **11,000 x g**. Discard the filtrate.
- b. Add **500 μ L Wash Buffer WB2** (with ethanol added). Centrifuge for **1 min** at **11,000 x g**. Discard the filtrate.
- c. Centrifuge for **2 min** at **11,000 x g**.

Note: Residual ethanol is removed during this step.

5. Elute DNA

- a. Place the Spin Column into a 1.5 mL microcentrifuge. Add **100 μ L Elution Buffer**. Incubate at room temperature for 1 min. Centrifuge **1 min** at **11,000 x g**.

Protocol for DNA Purification from Large Volumes of Whole Blood

1. Sample Preparation

- a. Add 1mL of ice cold nuclease free water to 500 μ L of whole blood, mix thoroughly by vortexing or pipetting.
- b. Incubate the sample for 5 min at room temperature.
- c. Centrifuge for 5 min at 3,000 rpm.
- d. Discard the supernatant.
- e. Resuspend the pellet in 200 μ L of 1 x PBS.
- f. Proceed to step 1 of the Whole Blood Genomic DNA Purification Protocol.
Note. For purification of DNA from samples exceeding the standard 200 μ L volume, it is necessary to burst red blood cells prior to performing the cell lysis step. Up to 500 μ L of mammalian blood can be processed using following protocol

Protocol for DNA Purification from Nucleated Blood

1. Sample Preparation

- a. Take 2-10 μ L of nucleated blood.
- b. Adjust the volume to 200 μ L with 1xPBS
- c. Proceed to step 1 of the Whole Blood Genomic DNA Purification Protocol
Note: Nucleated avian or fish blood contains very large amounts of genomic DNA and therefore the volume of the starting material has to be scaled down. The DNA purification procedure follows the same protocol as mammalian blood, except that 2-10 μ L of blood are used per purification.

Protocol DNA Purification from Buccal Swabs

1. Sample Preparation

- a. To collect a sample, scrape the swab 5-6 times against the inside cheek
- b. Swirl the swab for 30-60 s in 200 μ L of 1xPBS.
- c. Go to step 1 of the standard Whole Blood Genomic DNA Purification Protocol

Protocol for DNA Purification from Bone Marrow

1. Sample Preparation

- a. Harvest 25-200 μL of fresh or frozen bone marrow.
- b. Adjust the volume to 200 μL with 1 \times PBS
- c. Proceed to step 1 of the Whole Blood Genomic DNA Purification Protocol

Protocol for DNA Purification from Buffy Coat

1. Sample Preparation

- a. Centrifuge 1.5 mL of whole blood at 5,000 rpm for 10 min at room temperature. Three layers should be visible.
- b. Remove upper clear layer by aspiration
- c. Collect approximately 200 μL of intermediate layer using an automatic pipette.
Note: If necessary, adjust the volume to 200 μL with 1 \times PBS
- d. Proceed to step 1 of the Whole Blood Genomic DNA Purification Protocol.

Protocol for DNA Purification from Bone Marrow

1. Sample Preparation

- a. Add 0.5 mL of 0.5 M EDTA to 4.5 mL of urine (final concentration 50 mM).
- b. Centrifuge 10 min at 3000 rpm.
- c. Discard the supernatant.
- d. Resuspend the pellet in 200 μL 1 \times PBS.
- e. Proceed to step 1 of the Whole Blood Genomic DNA Purification Protocol

Protocol for DNA Purification from Dried Blood Spots

1. Sample Preparation

- a. Cut out the section of filter containing the dried blood sample and place into a microcentrifuge tube.
- b. Add 200 μL of 1 \times PBS and incubate 5-10 min at room temperature.
- c. Proceed to step 1 of the Whole Blood Genomic DNA Purification Protocol
Note: Taqgen[®] Tissue DNA kit ensures the efficient purification for the Dried Blood cards.